

# Product Datasheet

## Anti-Human CD49d Monoclonal Antibody, FITC Conjugated Clone PS/2 853.441.010

Product Name	Anti-Human CD49d Monoclonal Antibody, FITC Conjugated Clone PS/2
Catalog Number	853.441.010
Description	Mouse monoclonal
Species/Host	Rat
Conjugation	FITC
Tested Applications	FC
Brand	Diaclone
Form/Appearance	Phosphate buffered saline, 0,09% sodium azide, pH 7.2, containing 1% BSA
Storage	Store at +4°C protected from light.
Isotype	Rat IgG2b
Clonality	Monoclonal
Specificity	Reacts with the kD VLA-alpha-4 chain (lymphocyte homing receptor), which is associated with the CD29 molecule in the VLA-4 complex