

Product Datasheet

Anti-Human CD49d Monoclonal Antibody, Unconjugated Clone BU49 853.443.020

Product Name	Anti-Human CD49d Monoclonal Antibody, Unconjugated Clone BU49
Catalog Number	853.443.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IP, IHC
Immunogen	Hyatt myeloma cell line
Brand	Diaclone
Form/Appearance	Phosphate buffered saline, 0,09% sodium azide, pH 7.4, containing 1% BSA
Storage	Store at +2-8°C. For longer storage freeze aliquots.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgG1
Clonality	Monoclonal
Specificity	Reacts with the kD VLA-alpha-4 chain (lymphocyte homing receptor), which is associated with the CD29 molecule in the VLA-4 complex
Hybridoma	NSO myeloma cell line