

# Product Datasheet

## Anti-Human CD57 Monoclonal Antibody, Unconjugated Clone NC1 853.463.020

Product Name	Anti-Human CD57 Monoclonal Antibody, Unconjugated Clone NC1
Catalog Number	853.463.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC, IHC
Immunogen	Quail ciliary ganglia cells
Brand	Diaclone
Form/Appearance	20 mM TRIS and 150 mM NaCl with 1% BSA, pH 8.0. Preservative: Mertiolate (0.001%)
Storage	Store at +2-8°C. For longer storage freeze aliquots.
Note	We recommend that each laboratory determine an optimum working titre for use in its particular application
Isotype	IgM
Clonality	Monoclonal
Source	Quail ciliary ganglia cells
Specificity	Directed against the CD57-antigen, which is expressed on NK-cells and subset T lymphocytes
Hybridoma	Myeloma SP2/O Ag 1.4 x Balb/c spleen cells