

Product Datasheet

Anti-Human CD178 Monoclonal Antibody, Azide Free Clone B-B34 853.590.005

Product Name	Anti-Human CD178 Monoclonal Antibody, Azide Free Clone B-B34
Catalog Number	853.590.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	FasL transfected 293T cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises APO1L, Fas ligand (FAS L), a 40 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells