

## Product Datasheet

### Serum, Single Donor 991-24-PS

<b>Product Name</b>	Serum, Single Donor
<b>Catalog Number</b>	991-24-PS
<b>Description</b>	Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria.
<b>Tested Applications</b>	QC Materials
<b>Alternative Names</b>	Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Source</b>	Single Human Donor
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Preservatives</b>	None added, inquire if needed