

Product Datasheet

Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-C48 853.620.005

Product Name	Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-C48
Catalog Number	853.620.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA, ELISpot
Immunogen	Natural human Perforin
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes human Perforin, a 70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells