

## Product Datasheet

### Buccal Swab, Single Donor 991-27-S

<b>Product Name</b>	Buccal Swab, Single Donor
<b>Catalog Number</b>	991-27-S
<b>Description</b>	Each sample is equivalent to 2 swabs per donor. Collection performed a minimum of 30 minutes after the last eating/drinking event. Inquire for larger donor counts, swab counts, or custom criteria.
<b>Tested Applications</b>	QC Materials
<b>Alternative Names</b>	Cheek swab, oral swab, mouth swab, cheek, oral, mouth, salivary, buccal
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Swabs
<b>Storage</b>	-20°C
<b>Source</b>	Single Human Donor
<b>Shelf Life</b>	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
<b>Preservatives</b>	None added, inquire if needed