

Product Datasheet

Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-D30 853.630

Product Name	Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-D30
Catalog Number	853.630
Description	Mouse monoclonal antibody. Source: Myeloma X63/AG.8653 x Balb/c spleen cells
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	U266 cell line
Alternative Names	SYNDECAN-1
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	Recognises Syndecan-1, the heparan sulfate proteoglycan, a 65-70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells