

Product Datasheet

Perianal Swab, Single Donor, Male 991-34-M

Product Name	Perianal Swab, Single Donor, Male
Catalog Number	991-34-M
Description	Each sample is equivalent to 2 swabs per donor. Inquire for larger donor counts, swab counts, or custom criteria.
Tested Applications	QC Materials
Alternative Names	Perineum swabs, perineum, anorectal, perianal
Brand	Lee Biosolutions
Form/Appearance	Swabs
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed