

Product Datasheet

Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-D30 853.630.005

Product Name	Anti-Human CD138 Monoclonal Antibody, Azide Free Clone B-D30
Catalog Number	853.630.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	U266 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Syndecan-1, the heparan sulfate proteoglycan, a 65-70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells