

Product Datasheet

Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-E48 853.640.000

Product Name	Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-E48
Catalog Number	853.640.000
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA, ELISpot
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Source	Natural human Perforin
Specificity	Recognizes human Perforin, a 70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Natural human Perforin