

Product Datasheet

Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-E48 853.640.005

Product Name Anti-Human Perforin Monoclonal Antibody, Azide Free Clone B-E48

Catalog Number 853.640.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications ELISA, ELISpot

Immunogen Natural human Perforin

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Store at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognizes human Perforin, a 70 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells