

# Product Datasheet

## Plasma, Single Donor, Anti-Thyroglobulin (anti-TG) Positive 991-58-S-ATG

|                     |   |
|---------------------|---|
| Product Name        | Plasma, Single Donor, Anti-Thyroglobulin (anti-TG) Positive   |
| Catalog Number      | 991-58-S-ATG  |
| Tested Applications | QC Materials  |
| Alternative Names   | Blood plasma, plasma, whole plasma, fraction, plasma fraction, anticoagulated, thyroglobulin, anti-thyroglobulin, anti-TG, TG, thyroid  |
| Brand               | Lee Biosolutions  |
| Form/Appearance     | Liquid  |
| Concentration       | anti-TG: >= 1,000 IU/mL (Immunoassay)   |
| Storage             | -20°C   |
| Source              | Single Human Donor  |
| Shelf Life          | Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |
| Biological Activity | anti-TG: <u>≥</u> 1,000 IU/mL (Immunoassay)   |
| Preservatives       | Typically EDTA (anticoagulant), inquire for other options   |