

Product Datasheet

Plasma, Single Donor, Human Anti-Mouse Antibody (HAMA) Positive 991-58-S-HAMA

Product Name	Plasma, Single Donor, Human Anti-Mouse Antibody (HAMA) Positive
Catalog Number	991-58-S-HAMA
Tested Applications	QC Materials
Alternative Names	Blood plasma, plasma, whole plasma, fraction, plasma fraction, anticoagulated, HAMA, mouse, anti-mouse
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	HAMA: Typically ≥ 100 ng/mL (Immunoassay)
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	Typically EDTA (anticoagulant), inquire for other options