

Product Datasheet

Anti-h D-Dimer 1403 SPTN-5 100228

| | |
|----------------------------|--|
| Product Name | Anti-h D-Dimer 1403 SPTN-5 |
| Catalog Number | 100228 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | FIA, IT |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8°C |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Epitope | N/D |
| Purity | ≥ 95 % |
| Associated Products | Native D-Dimer, Lee Biosolutions, Cat. 200-09, 200-12 and 200-13 |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.7-7.6 |
| Cross Reactivity | Recognizes human fibrinogen |
| Specificity | Antibody recognizes human D-dimer |
| Shelf Life | 36 months |