

Product Datasheet

Anti-h D-Dimer 1403 SPTN-5 100228

Product Name	Anti-h D-Dimer 1403 SPTN-5
Catalog Number	100228
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Species/Host	Mouse
Tested Applications	FIA, IT
Alternative Names	D dimer
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Epitope	Not determined
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.7-7.6
Cross Reactivity	Recognizes human fibrinogen
Specificity	Antibody recognizes human D-dimer
Shelf Life	36 months