

Product Datasheet

Anti-Human p53 Monoclonal Antibody, Azide Free Clone B-F34 853.660.005

Product Name	Anti-Human p53 Monoclonal Antibody, Azide Free Clone B-F34
Catalog Number	853.660.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant p53 protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises wild type and mutant p53, the epitope is within residues 18-30 of human p53
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells