

Product Datasheet

Luteinizing Hormone (LH), Human Pituitary Gland 996-31

| | |
|----------------------------|---|
| Product Name | Luteinizing Hormone (LH), Human Pituitary Gland |
| Catalog Number | 996-31 |
| Alternative Names | LH, Lutropin, Interstitial Cell-Stimulating Hormone |
| Brand | Lee Biosolutions |
| Form/Appearance | Lyophilized |
| Purity | ≥ 98% (SDS-PAGE) |
| Source | Human Pituitary Glands |
| Molecular Weight | ~30,000 |
| Shelf Life | 5 years |
| Biological Activity | ≥ 8,500 IU/mg |
| Content | FSH: ≤ 18.75 IU/mg (All assayed) ≤ 0.1% w/w) TSH: ≤ 0.0425 IU/mg PRL: ≤ 0.005 mg/mg hGH: ≤ 0.005 mg/mg |