

Product Datasheet

Anti-Human IFN γ Monoclonal Antibody, Azide Free Clone B-G1 853.670.005

Product Name	Anti-Human IFN γ Monoclonal Antibody, Azide Free Clone B-G1
Catalog Number	853.670.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA, ELISpot
Immunogen	Recombinant human IFN-g
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 μ m. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IFN-g
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells