

Product Datasheet

Immunoglobulin A (IgA), anti-Human Polyclonal Antibody GA-80A

Product Name	Immunoglobulin A (IgA), anti-Human Polyclonal Antibody
Catalog Number	GA-80A
Species/Host	Goat
Immunogen	High purity native human plasma Immunoglobulin A
Alternative Names	anti-IgA
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1 - 2 mg/mL
Storage	+2-8 °C
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity Purified
Shelf Life	2 years