

Product Datasheet

Albumin (ALB), Goat Polyclonal anti-Human Antibody - Affinity Purified GAL-80A

Product Name	Albumin (ALB), Goat Polyclonal anti-Human Antibody - Affinity Purified
Catalog Number	GAL-80A
Species/Host	Goat
Immunogen	High purity native human serum Albumin
Alternative Names	Albumin pAb, HSA Ab, ALB Ab
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	+2-8 °C
Clonality	Polyclonal
Purity	Affinity Purified
Shelf Life	2 years