

Product Datasheet

Complement Component 3 (C3), Goat Polyclonal Antibody GC3-80A

Product Name	Complement Component 3 (C3), Goat Polyclonal Antibody
Catalog Number	GC3-80A
Species/Host	Goat
Immunogen	High purity human plasma Complement Component 3
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear solution
Storage	2-8° C
Isotype	IgG
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years