

## Product Datasheet

### C4 Polyclonal Goat Anti-Human Affinity Purified Antibody GC4-80A

<b>Product Name</b>	C4 Polyclonal Goat Anti-Human Affinity Purified Antibody
<b>Catalog Number</b>	GC4-80A
<b>Description</b>	Purified Polyclonal Antibody
<b>Species/Host</b>	Goat
<b>Immunogen</b>	High purity native human plasma Complement 4 (C4)
<b>Alternative Names</b>	B-type Natriuretic Peptide
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity Purified
<b>Shelf Life</b>	2 years