

Product Datasheet

Polyclonal anti-Human Ceruloplasmin Antibody GC-80A

Product Name	Polyclonal anti-Human Ceruloplasmin Antibody
Catalog Number	GC-80A
Species/Host	Goat
Tested Applications	ELISA, Blotting, IP
Immunogen	High purity, native human Ceruloplasmin
Alternative Names	anti-Cp, anti-Ceruloplasmin
Brand	Lee Biosolutions
Form/Appearance	Liquid Clear solution
Storage	2-8° C
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years