

Product Datasheet

Creatine Kinase MM (CK-MM), anti-Human Polyclonal Antibody GCKM-80A

Product Name	Creatine Kinase MM (CK-MM), anti-Human Polyclonal Antibody
Catalog Number	GCKM-80A
Species/Host	Goat
Immunogen	High purity native human heart Creatine Kinase MM
Alternative Names	CKMM, CK-MM
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8° C
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years