

## Product Datasheet

### C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified GCRP-80A

<b>Product Name</b>	C-Reactive Protein (CRP), Polyclonal Antibody Anti-Human Affinity Purified
<b>Catalog Number</b>	GCRP-80A
<b>Species/Host</b>	Goat
<b>Immunogen</b>	Highly purified human C-Reactive Protein (CRP)
<b>Brand</b>	Lee Biosolutions
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 - 2 mg/mL
<b>Storage</b>	+2-8 °C
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity Purified
<b>Shelf Life</b>	2 years