

Product Datasheet

Goat Cystatin C Polyclonal Affinity Purified GCYS-80A

Product Name	Goat Cystatin C Polyclonal Affinity Purified
Catalog Number	GCYS-80A
Species/Host	Goat
Immunogen	High purity human Cystatin C
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1 - 2 mg/mL (A_{280} , $E_{1\text{mg/mL}}^{1\text{cm}}_{280\text{nm}} = 1.4$)
Storage	2-8° C
Clonality	Polyclonal
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 Years