

Product Datasheet

Neuron Specific Enolase (NSE), Polyclonal Antibody GNEN-80A-Z

Product Name	Neuron Specific Enolase (NSE), Polyclonal Antibody
Catalog Number	GNEN-80A-Z
Species/Host	Goat
Immunogen	High purity human neuron specific enolase
Alternative Names	anti-NSE
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8° C
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years