

Product Datasheet

NGAL, Goat anti-Human Polyclonal Antibody GNGAL-80A-Z

Product Name	NGAL, Goat anti-Human Polyclonal Antibody
Catalog Number	GNGAL-80A-Z
Species/Host	Goat
Immunogen	High purity recombinant human NGAL Antigen
Alternative Names	Lipocalin-2, LCN2, oncogene 24p3, neutrophil gelatinase-associated lipocalin
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1 - 2 mg/mL
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	1 year