

Product Datasheet

Prealbumin (Transthyretin, TTR), anti-Human Polyclonal Antib GPRE-80A

Product Name	Prealbumin (Transthyretin, TTR), anti-Human Polyclonal Antib
Catalog Number	GPRE-80A
Species/Host	Goat
Immunogen	Highly purified human prealbumin
Alternative Names	Transthyretin, Thyroxine Binding Prealbumin, ATTR, TBPA
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	1 - 2 mg/mL (A _{280nm} , E _{0.1%} = 1.4)
Storage	2-8° C
Isotype	IgG
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years