

Product Datasheet

Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat GRBP-80A

Product Name	Retinol Binding Protein (RBP), Polyclonal Anti-Human Affinity Purified Host Goat
Catalog Number	GRBP-80A
Species/Host	Goat
Immunogen	High purity, native human RBP
Alternative Names	anti-RBP, anti-Retinol Binding Protein
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8° C
Isotype	IgG
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years