

Product Datasheet

Ferritin, Polyclonal Antibody Anti-Human Affinity Purified Purified RF-80A

Product Name	Ferritin, Polyclonal Antibody Anti-Human Affinity Purified
Catalog Number	RF-80A
Description	Affinity Purified Polyclonal Antibody
Species/Host	Rabbit
Immunogen	High purity native human liver Ferritin
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	+2-8 °C
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years