

Product Datasheet

Anti-Human L1CAM Monoclonal Antibody, Azide Free Clone B-R48 857.370.005

Product Name	Anti-Human L1CAM Monoclonal Antibody, Azide Free Clone B-R48
Catalog Number	857.370.005
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant human L1CAM/Fc
Alternative Names	CD171 / NCAM-L1
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Clone: B-R48
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human L1CAM
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells