

Product Datasheet

Anti-SARS-CoV-2 RBD domain Monoclonal Antibody, Azide Free Clone B-K45 857.380.005

Product Name	Anti-SARS-CoV-2 RBD domain Monoclonal Antibody, Azide Free Clone B-K45
Catalog Number	857.380.005
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 RBD domain
Alternative Names	Receptor binding domain Spike protein / RBD S1 protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free.
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Clone: B-K45
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	RBD domain from Spike protein (aa319-541)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells