

## Product Datasheet

### Anti-Human LRRC15 Monoclonal Antibody, Azide Free Clone B-G53 857.400

<b>Product Name</b>	Anti-Human LRRC15 Monoclonal Antibody, Azide Free Clone B-G53
<b>Catalog Number</b>	857.400
<b>Description</b>	Source: Myeloma X63/AG.8653 x Balb/c lymph nodes cells
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Azide free
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human LRRC15
<b>Alternative Names</b>	Leucine rich repeat containing 15   LIB protein   hLib
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/ml (+/- 10 %)
<b>Storage</b>	Stable at +2-8 °C for 12 months. DO NOT FREEZE
<b>Note</b>	Clone: B-G53
<b>Isotype</b>	Mouse IgG2b Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
<b>Specificity</b>	Recognizes natural and recombinant human LRRC15
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c node cells