

Product Datasheet

Human IFN gamma EasySplit ELISpot Kit 856.051.005PCS

Product Name	Human IFN gamma EasySplit ELISpot Kit
Catalog Number	856.051.005PCS
Description	Diaclone EasySplit ELISpot kits include Pre-coated strip plate, Detection biotinylated antibody, Streptavidin Alkaline Phosphatase Conjugate, Bovine Serum Albumin, Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
Specificity	Recognizes natural human IFN-g