

Product Datasheet

Murine IFN gamma EasySplit ELISpot Kit 862.031.002PCS

Product Name	Murine IFN gamma EasySplit ELISpot Kit
Catalog Number	862.031.002PCS
Description	Diaclone EasySplit ELISpot kits include Pre-coated strip plate, Detection biotinylated antibody, Streptavidin Alkaline Phosphatase Conjugate, Bovine Serum Albumin, Ready-to-use BCIP/NBT Substrate Buffer
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Cross Reactivity	No cross reactivity with other murine cytokines
Specificity	Recognizes natural murine IFN-g