

Product Datasheet

Human RANTES ELISA Kit 873.110.048

Product Name	Human RANTES ELISA Kit
Catalog Number	873.110.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	250 – 6.25 pg/ml
Sensitivity	1.2 pg/ml
Tested Applications	ELISA
Alternative Names	CCL5, C-C motif chemokine ligand 5
Brand	Diaclone
Cross Reactivity	No cross reactivity with other cytokines tested
Specificity	Recognizes both natural and recombinant RANTES
Sample Types	Cell culture supernatants, serum and plasma