

# Product Datasheet

## Human ST2 ELISA Kit 873.120.048

Product Name	Human ST2 ELISA Kit
Catalog Number	873.120.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	5000 – 156.25 pg/ml
Sensitivity	38 pg/ml
Tested Applications	ELISA
Alternative Names	IL-1R4,IL-33 receptor,IL-1RL1
Brand	Diaclone
Cross Reactivity	No cross reactivity with other cytokines tested
Specificity	Recognizes both natural and recombinant ST2
Sample Types	Cell culture supernatants, serum and plasma