

Product Datasheet

Human EpCAM ELISA kit 873.130

Product Name	Human EpCAM ELISA kit
Catalog Number	873.130
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	20 - 0,625 ng/ml
Sensitivity	0.27 ng/ml
Alternative Names	CD326, TROP1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other cytokines tested
Specificity	Recognizes both natural and recombinant EpCAM
Sample Types	Cell culture supernatants, serum and plasma