

Product Datasheet

Human EpCam ELISA kit 873.130.192

Product Name	Human EpCam ELISA kit
Catalog Number	873.130.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	20 – 0,625 ng/ml
Sensitivity	0.27 ng/ml
Tested Applications	ELISA
Alternative Names	CD326, TROP1
Brand	Diaclone
Cross Reactivity	No cross reactivity with other cytokines tested
Specificity	Recognizes both natural and recombinant EpCAM
Sample Types	Cell culture supernatants, serum and plasma