

Product Datasheet

Recombinant Human B Cell Maturation Protein (BCMA) 715-H29-100

Product Name	Recombinant Human B Cell Maturation Protein (BCMA)
Catalog Number	715-H29-100
Species/Host	HEK cells
Alternative Names	CD269, TNFRSF17, BCMA tumor necrosis factor receptor superfamily member 17
Brand	Diaclone
Form/Appearance	Liquid. Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm
Storage	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
Purity	>90% by SDS PAGE gel
Sequence	ECD Leu2-Ala54
Tag	Silent Human IgG1 Fc (Pro100-Lys330)