

Product Datasheet

Anti-Human BCMA Monoclonal Antibody, Azide Free Clone B-B54 857.410.005

Product Name	Anti-Human BCMA Monoclonal Antibody, Azide Free Clone B-B54
Catalog Number	857.410.005
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC
Immunogen	Recombinant human BCMA ECD with Fc Human IgG1 domain
Alternative Names	CD269 / TNFRSF17 / B cell maturation antigen
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Clone: B-B54
Isotype	Mouse IgG1 Kappa light chain
Source	Recombinant antibody from a phage display library
Specificity	Recognizes both natural and recombinant human BCMA