

# Product Datasheet

## Anti-Human BCMA Monoclonal Antibody, Unconjugated Clone B-B54 857.413.020

Product Name	Anti-Human BCMA Monoclonal Antibody, Unconjugated Clone B-B54
Catalog Number	857.413.020
Description	Source: Recombinant antibody from a phage display library
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human BCMA ECD with Fc Human IgG1 domain
Alternative Names	CD269   TNFRSF17   B cell maturation antigen
Brand	Diaclone
Form/Appearance	Liquid
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Note	Clone: B-B54
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Recombinant antibody from a phage display library
Buffer	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Specificity	Recognizes both natural and recombinant human BCMA