

Product Datasheet

Anti-Human CX3CL1 Monoclonal Antibody, Azide Free Clone B-F51 857.420.005

Product Name	Anti-Human CX3CL1 Monoclonal Antibody, Azide Free Clone B-F51
Catalog Number	857.420.005
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant human CX3CL1
Alternative Names	Fractalkine / C-X3-C motif chemokine ligand 1
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE
Note	Clone: B-F51
Isotype	Mouse IgG2a Kappa light chain
Specificity	Recognize natural and recombinant human CX3CL1
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells