

## Product Datasheet

### Cotinine-4-BSA Antigen 170047

|                        |   |
|------------------------|---|
| <b>Product Name</b>    | Cotinine-4-BSA Antigen  |
| <b>Catalog Number</b>  | 170047  |
| <b>Description</b>     | Synthetic Cotinine, BSA conjugated  |
| <b>Brand</b>           | Medix Biochemica  |
| <b>Form/Appearance</b> | Liquid; Antigen-BSA conjugate   |
| <b>Storage</b>         | Short term storage ( $\leq 2$ weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.                         |
| <b>Note</b>            | This product was previously available with catalog code P82-99-04A. It has been produced using the same conjugation protocol. |
| <b>Buffer</b>          | Distilled water   |
| <b>Matched Pairs</b>   | PP194IF or PP198IF (detection)  |
| <b>Shelf Life</b>      | 60 months   |