

Product Datasheet

Cotinine-4-BSA Antigen 170047

Product Name	Cotinine-4-BSA Antigen
Catalog Number	170047
Description	Synthetic Cotinine, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short term storage (≤ 2 weeks): 2-8 °C. Long term:-20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with catalog code P82-99-04A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	PP194IF or PP198IF (detection)
Shelf Life	60 months