

Product Datasheet

Clonazepam-BSA Antigen 170064

Product Name	Clonazepam-BSA Antigen
Catalog Number	170064
Description	Synthetic Clonazepam, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog codes P83-99-21A and LA270. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	100912 (detection)
Shelf Life	60 months