

Product Datasheet

Seurm, Single Donor, Prolactin (PRL) 991-24-PS-PRL

Product Name	Seurm, Single Donor, Prolactin (PRL)
Catalog Number	991-24-PS-PRL
Description	Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria.
Tested Applications	QC Materials
Alternative Names	Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated, prolactin, lactotropin, mammotropin, PRL
Brand	Lee Biosolutions
Form/Appearance	Liquid
Concentration	Reported in ng/mL
Storage	-20°C
Source	Single Human Donor
Shelf Life	Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time.
Preservatives	None added, inquire if needed