

Product Datasheet

Serum, Single Donor, Sex Hormone Binding Globulin (SHBG) 991-24-PS-SHBG

| | |
|---------------------|---|
| Product Name | Serum, Single Donor, Sex Hormone Binding Globulin (SHBG) |
| Catalog Number | 991-24-PS-SHBG |
| Description | Each sample is a minimum of 1mL per donor. Inquire for larger volumes per donor or custom criteria. |
| Tested Applications | QC Materials |
| Alternative Names | Serum, fraction, clotted, post-coagulation, coagulated, powder, sera, blood serum, defibrinated, SHBG, sex hormone, sex steroid, SSBG |
| Brand | Lee Biosolutions |
| Form/Appearance | Liquid |
| Concentration | Reported in nmol/L |
| Storage | -20°C |
| Source | Single Human Donor |
| Shelf Life | Due to differing stabilities of various biomarkers in this complex biological material, it is not possible to assign a single expiration date that is relevant to all biomarkers. The end user should monitor the biomarkers of interest to verify suitability over time. |
| Preservatives | None added, inquire if needed |